Aquafluor®

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Page 3 - en: Instructions for use

Page 4 - fr: Manuel d'instruction

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Strana 7 - bs: Uputstvo za upotrebu

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Composition

Qualitative data:

Sodium fluoride, bis-(hydroxyethyl)-aminopropyl-N-hydroxyethyl-octadecylamine dihydrofluoride, aqueous methacrylate dispersion

Quantitative data:

Sodium fluoride, bis-(hydroxyethyl)-aminopropyl-N-hydroxyethyl-octadecylamine dihydrofluoride 0.15 % calculated as pure fluoride

Intended purpose of the product

Fluoridation of the tooth enamel.

Indications

Aquafluor is a varnish for fluoridation of the tooth enamel and desensitisation of sensitive teeth.

Contraindications

Do not use in patients with a known allergy to methacrylates.

Performance characteristics of the product

After varnish formation (application of dispersion and dry blowing) on the enamel, the varnish softens and releases fluoride.

Specific precautions

None

Properties

Aquafluor is an air-drying dispersion that forms a well-adhering film on enamel and dentin. The adherent film softens when it comes into contact with saliva or water for a longer time. In the process, the fluorides contained in the polymer matrix can diffuse to the dentin or enamel surface and exert their effect there. Due to this behaviour, minimal amounts of fluoride are completely sufficient for optimal fluoridation. Once the varnish is completely soaked, the fluoride release is complete. The varnish matrix then detaches from the enamel or dentin surface.

Application

Clean surfaces to be treated and, if possible, dry them. Transfer Aquafluor drop by drop to the brush (it is not necessary to shake up the aqueous dispersion before removal) and touch surfaces thinly. Dry the dispersion carefully and extensively, preferably with warm air, until the polymer film has completely formed and is dry. Do not brush teeth for about 24 hours after treatment. When used correctly, the polymer layer will not come off until it is completely soaked and maximum fluoride release has occurred.

NOTE

The application aids included with the product (cannulas, disposable brushes, etc.) are only for single use on a single patient!

Duration of application

For caries prophylaxis, treatment should be carried out 2–3 times a year (or 3–4 times if there is an increased susceptibility to caries).

Warnings**Residual risks**

The processing instructions and precautionary measures described in the instructions for use of the product must be observed in every instance. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, please inform us in every case, describing the accompanying circumstances as precisely as possible.

Side effects

Undesirable side effects of this medical product are extremely rare, provided the product is processed and used properly. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible. We follow up on every indication.

No risks in the event of an overdose are known of.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

The Product should be stored between 0° C - 22°C.

Minimum shelf life

Do not use after the expiry date.

Variants

REF WP4005 1 x bottle with 5.0 ml liquid + 50 brushes

REF WP4004 1 x bottle with 10.0 ml liquid

REF WP4949 1 x bottle with 5.0 ml liquid, Strawberry + 50 brushes

REF WP4948 1 x bottle with 10.0 ml liquid, Strawberry

Disposal

Dilute unused liquid with water and dispose of with household waste. Cured residues can also be disposed of with household waste.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

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