Competence Flow®

- Seite 2 - de: Gebrauchsanweisung
Page 3 - en: Instructions for use
Page 4 - fr: Manuel d'instruction
Página 5 - es: Instrucciones de uso
Pagina 6 - it: Istruzioni per l'uso
Strana 7 - bs: Uputstvo za upotrebu
Pagină 8 - ro: Instrucțiuni de utilizare
Sayfa 9 - tr: Kullanma kılavuzu
Strona 10 - pl: Instrukcja użytkowania
Página 11 - pt: Manual de instruções
Стр. 12 - ru: Инструкция по применению

We know how

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Patientenzielgruppe

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Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended user

Dentist

Composition

Qualitative composition:

1,4-butanediol dimethacrylate, bis-GMA, aliphatic urethane dimethacrylate, butylhydroxytoluene, camphorquinone, silicon dioxide, glass ceramic, iron oxide pigments, 2-(2H-benzotriazol-2yl)-p-cresol, 2-ethylhexyl 4-(dimethylamino)benzoate, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

Quantitative composition:

1,4-butanediol dimethacrylate 5%, aliphatic urethane dimethacrylate 22%, bis-GMA 18%, silicon dioxide 6%, activators and UV absorbers <1%, glass ceramic 48%, pigments <1%

Intended purpose of the product

Competence Flow is a light-curing, radiopaque, flowable composite filling material for the dental restoration of Class III, IV and V cavities in children, adolescents, adults and elderly patients by a dentist.

Indications

For fillings of cavity classes III, IV and V.

Contraindications

Do not use on patients with known allergies to the ingredients.
Do not use directly on the opened pulp.
Use is contraindicated if draining or the prescribed application technique is not possible.

Performance characteristics of the product

Competence Flow is a flowable composite that cures with light.

Specific precautions

None

Application

- Clean and dry the preparation. Determine the colour.
- Prepare the cavity according to general adhesive technique rules. Bevel the enamel edges.
- Etch the enamel margins with etching gel for 20 to 35 seconds.
- Rinse off etching gel completely.
- Dry with oil-free air. Avoid overdrying the dentine.
- Apply adhesive according to the instructions for use.
- Select the syringe in the desired colour. Remove syringe cap. Screw on the Luer-Lock needle.
- Competence Flow can be applied directly by pressing the syringe. To ensure complete hardening, apply Competence Flow in layers no thicker than 2 mm.
- Light-cure each layer for 30 to 40 seconds.
The halogen or LED light curing device should emit at 450 nm with an illumination intensity $\geq 600 \text{ mW/cm}^2$ and be checked regularly. Position the light as close as possible to the filling material. Do not remove the oxygen-inhibited layer, as it is required for bonding with the next layer. The use of a mould is recommended.
- Repeat steps 8 and 9 until the cavity is filled.
- Excess removal, contouring and finishing are carried out using diamond finishers and flexible discs with decreasing grain size.

Precautions

The user and patient should wear protective goggles before application, as well as when applying etching gel.

Avoid eye contact.

In case of eye contact, rinse immediately with plenty of water and seek medical advice.

Avoid contact with the skin.

In the event of skin contact, remove with a cloth and alcohol and then wash thoroughly with soap and water. In case of hypersensitivity of the skin or skin rash, discontinue use and consult a doctor.

The material contains titanium dioxide and glass ceramic, and it should only be ground outside the mouth under an intact fume hood.

Details

Delivery only to dentists and dental technicians or on their behalf. Keep out of the reach of children! Only dispose of completely emptied packaging for recycling.

Side effects

No undesirable side effects are known.

Interactions

Avoid materials containing eugenol, moisture and oily air, as these can prevent polymerisation at the point of contact.

Discolouration may occur if cationic mouthwashes, plaque disclosing agents or chlorhexidine are used.

Storage

Do not expose to direct sunlight or other heat sources. Store in a dry place.

The product should be stored between 0°C to 22°C.

Minimum shelf life

3 years. Do not use after the expiry date.

Variants

-  WP4197 3× syringes (A2, A3, A3.5), 1.0 ml + 9 needles
-  WP4199 1× syringe A1 with 1.0 ml + needles
-  WP4190 1× syringe A2 with 1.0 ml + needles
-  WP4191 1× syringe A3 with 1.0 ml + needles
-  WP4192 1× syringe A3.5 with 1.0 ml + needles
-  WP4350 1× syringe B1 with 1.0 ml + needles
-  WP4196 1× syringe B2 with 1.0 ml + needles
-  WP4398 1× syringe C2 with 1.0 ml + needles
-  WP4370 1× syringe A1 with 3.5 g + needles
-  WP4380 1× syringe A1 with 3.5 g + needles
-  WP4371 1× syringe A2 with 3.5 g + needles
-  WP4381 1× syringe A2 with 3.5 g + needles
-  WP4372 1× syringe A3 with 3.5 g + needles
-  WP4382 1× syringe A3 with 3.5 g + needles
-  WP4373 1× syringe A3.5 with 3.5 g + needles
-  WP4383 1× syringe A3.5 with 3.5 g + needles
-  WP4384 1× syringe B1 with 3.5 g + needles
-  WP4385 1× syringe B2 with 3.5 g + needles

Disposal

Any unused material can be disposed of in household waste after curing.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

We know how

dent a pharm
Produktionsgesellschaft mbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9225-0
Fax: +(49) 4123/9225-48
www.dent-a-pharm.de
info@dent-a-pharm.de

Willmann & Pein GmbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9228-0
Fax: +(49) 4123/9228-49
www.wp-dental.de
info@wp-dental.de

Made in Germany

H317 May cause allergic skin reactions
H318 Causes serious eye damage

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