Doublident®

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Page 3 - en: Instructions for use

Page 4 - fr: Manuel d'instruction

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Strana 7 - bs: Uputstvo za upotrebu

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Instructions for use

Patient target group

None, as Doublident does not come into contact with the human body.

Intended users

Dental technicians

Composition

Addition-curing silicones (vinyl polysiloxanes).

Intended purpose of the product

Silicone for model moulding.

Indications

Doublident is a moulding silicone for use in dental laboratories. The base/catalyst ratio is 1:1, which makes batching and mixing easier. The duplicating silicone can be mixed manually or with duplicating equipment. It is very fluid and does not require a vacuum. Doublident is very accurate in capturing details; the dimensions are durable and stable. Doublident is contact-resistant to various types of resins, coatings and gypsums.

Application

Model preparation:

Before using Doublident, ensure that the model to be moulded is perfectly clean. When milling, carefully remove any greasy residue. Degrease with a solvent such as acetone. Wait for steam to evaporate, if used, and ensure that the model is completely dry. Select an appropriately sized duplication flask and secure the model in the duplication flask with wax.

Batching and mixing either manually and/or using vacuum technique:

- 1) Shake the bottles of the individual components before use. Measure and weigh equal amounts of catalyst (white) and base (coloured) in a graduated container,
 - 2) Mix the two components well until a uniform colour is achieved.
- Attention: Make sure that no unmixed residue or paint splashes remain on the container walls.

How to use:

Pour the Doublident slowly into the duplication flask from a height of approx. 30 cm if possible. If the quantity measured out is not sufficient to complete the duplication, wait until the silicone has completely set and then add more material in the same way. The material bonds perfectly without compromising the end result. Use a jet of compressed air between the model and the mould created. Do not apply force with a spatula or other objects and do not remove from the duplication flask or you might deform the mould. Doublident is compatible with all types of gypsums, veneers or casting resins. We recommend using Doublident at room temperature; low temperatures lead to a longer setting time, higher temperatures shorten the setting time. Please refer to the technical data table for mixing, processing and setting times.

Important recommendations:

You must maintain the correct ratios of 1:1 in order to obtain the correct times and not to affect the final characteristics of the product. If using materials of different kinds for the die-spacers, resins, etc., you must wait for these products to stabilise completely to prevent a risk of inhibition of the polymerisation reaction. The surfaces with which the material comes into contact must be clean, grease-free and dry. Close the bottles after use.

Warnings**Residual risks**

None

Side effects

No undesirable side effects are known.
There are no known risks of overdose.

Interactions

There is no known interaction of the product with other medical products, pharmaceuticals or other substances.

Storage

Guaranteed for 24 months if the material is stored between 5° and 27°C (41°-80°F).

Minimum shelf life

Do not use after the expiry date.

Variants

[REF]WP5080 – Doublident orange 2 x 1000 ml Cat & Base

Disposal

Dispose of in accordance with legal regulations.

Any serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user is based.

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Willmann & Pein GmbH