D

Ideal Putty®

Gebrauchsinformation

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Ideal Putty®

Instructions for use

Ideal Putty is used as a base impression material for the (two step) double impression technique in conjunction with Ideal Light. Ideal Putty is kneadable precision impression material based on addition curing vinyl silicone.

Ideal Putty is slightly hydrophilic and compatible with moist surfaces. Ideal Putty is easy to mix and not sticky.

Field of Application

Perforated or non-perforated metal or plastic trays may be used with Ideal Putty. In order to improve adhesion of Ideal Putty on the tray, an adhesive could be applied.

Use equal parts (volume or weight) of base paste and catalyst paste. Take the base paste with the coloured spoon and the catalyst paste with the white spoon.

Knead until the colour is homogeneous (at least 45 seconds). The ready kneaded material is loaded into the prepared tray.

Insert the loaded tray in the mouth not later than 2 minutes. Recommended time in mouth: 2 min. The impression should be stored at normal room temperature (max. 25°C/ 77°F).

Note: Times mentioned for 23°C. Higher temperature shorten, lower temperature prolongs these times. The impression should be stored at normal room temperature (max. 25°C).

Casting

The impression can be worked out 60 minutes after removal from the mouth.

It should be stored at room-temperature (max. 77°F).

Electroplating

All impressions can be copper- or silver-plated.

Disinfection

Impressions can be disinfected with standard disinfection solutions.

Side-effects

Side-effects with this medical product are extremely rare if it is properly prepared and applied. Should you encounter undesirable side-effects - even in case of doubt - then please contact us and give a precise description of the circumstances and symptoms.

Processing instructions

Ensure that the processing instructions and precautionary measures described in the supplied product information are always strictly observed.

Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the circumstances.

Storage

Shelf life at least 2 years, when stored at 23°C.

Do not store above 25°C.

Additional notes

Do not interchange lids of pots.

Do not interchange spoons.

Do not use latex gloves!

Recommended are gloves based on polyethylene.

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Willmann & Pein GmbH