TR

## Aquaflur®

### Kullanma talimatı

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## Aquaflur®

### Directions for Use

Aquaflur is a varnish for the fluoridation of teeth.

#### Form of administration

A solvent-free fluoride varnish on an aqueous acrylate resin basis for brushed application.

#### Active ingredients

Sodium fluoride, bis-(hydroxyethyl)-aminopropyl-N-hydroxyethyl-oktadecylamine-dihydrofluoride totalling 0,15% calculated as pure fluoride.

#### Properties

Aquaflur is an air-drying emulsion that forms a film that adheres firmly to the tooth enamel and dentine. The adhering film soaks up moisture when in prolonged contact with saliva or water. The fluorides contained in the polymer matrix diffuse to the dentine or enamel surface where they develop their action. As a result of these properties a minimum amount of fluoride is quite sufficient to achieve optimal fluoridation. Fluoride release ends when the enamel has become completely soaked through after which the varnish matrix becomes separated from the enamel or dentine surface.

#### Manner of application

Clean and dry the surface of the teeth. Drip Aquaflur on a brush and brush thinly on the surface of the teeth. The aqueous emulsion does not have to be shaken prior to application. Dry the emulsion with plenty of air until the polymer film has fully set and is dry. After application do not brush the teeth for approx. 24 hours. If correctly applied, the polymer layer will only become detached after it has become completely soaked through for maximum fluoride release.

#### Side-effects

Undesirable side-effects of the

medical product are extremely rare if it is properly applied. In the event that undesirable side-effects do appear - even if in case of doubt - please inform us and give a precise description of the accompanying circumstances and symptoms. We will definitely follow up all notifications.

#### Duration of application

Treat 2-3 times per year for caries prevention, and 3-4 times if there is increased susceptibility to caries.

#### Application remarks

The application instructions and other precautionary measures described in the directions for the use of this product must be observed. If patients develop undesirable reactions or other effects, even though these instructions have been carefully observed, please notify us in all cases and give a precise description of the accompanying circumstances.

#### NOTE

Applicator-tools (tips, brushes) are only determinate for single-use at only one patient.

#### Package

WP4005 Bottle containing 5 ml fluid + 50 brushes  
WP4004 Bottle containing 10 ml

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