BH

## Competence Flow®

### Uputstvo za upotrebu

#### Opis

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Class III, IV and V restorations.

#### Recommended Procedures

- Clean out and dry preparation. Choose shade.
- Prepare the cavity according to the general rules for adhesive dentistry. Bevel enamel margins. Protect pulp in deep cavities, according to general rules of preparative dentistry.
- Etch enamel margins with Etch for 20 to 35 seconds.
- Rinse off etching gel completely.
- Dry with oil-free air. Avoid dentin desiccation.
- Apply the adhesive (Bonding) according to its instructions.
- Select the syringe containing the desired shade of Competence Flow. Remove syringe cap and insert an angled luer lock syringe tip. Lock tip into place.
- Pressing the syringe plunger dispenses Competence Flow.
- To ensure a complete cure do not place Competence Flow in layers exceeding 2mm.
- Light cure each layer of Competence Flow for 30 to 40 seconds.
- Repeat steps 8 and 9, until the cavity is filled.
- Remove excess, shape and finish with diamond burs and flexible disks of increasing grits.

#### Essential notes

Light curing units should have 450nm output and should be checked regularly. Place the light as close as possible to the restorative material. Do not disturb the oxygen inhibited layer because it is needed for bonding to the next layer. The use of a matrix is advised. Follow the manufacturer's instructions for use of any materials which are used with Competence Flow.

#### Safety Precautions

Placement is contraindicated, if a dry working area or the described application technique is not possible. Do not use on the exposed pulp or if an allergy to any component is known. Side effects are not known to date. In singular cases hypersensitivity to components of the material cannot be excluded. Avoid contact of paste with skin. If accidental contact occurs, wash immediately with soap and water. In case of contact with eyes rinse with copious amounts of water and consult a physician immediately. Avoid

materials containing eugenol, moisture and oily air; they can inhibit the polymerisation at the contact area. The use of cationic mouth washes, plaque indicators and chlorhexidine may lead to discolorations.

#### Side effects

No side effects are known to date. Hypersensitivity against components of the material cannot be excluded in every case. Avoid skin contact with the paste. In the event of inadvertent contact, wash immediately with soap and water. In the event of eye contact, immediately rinse with plenty of water and seek medical treatment. If you experience any undesired effects, even in case of doubt, we ask for notification in every case, describing the concomitants and syptoms as detailed as possible.

#### Processing instructions

Observe the processing instructions in the instructions for use included with the product and the precautionary measures under all circumstances. Should the patient still show an undesired reaction or if there are undesired effects, we ask for notification in every case,describing the concomitants as detailed as possible.

#### Compositon

Dental glass in an optimized BIS-GMA (oligomere Dimethacrylates) resin matrix, containing 63% inorganic fillers.

#### Storage and shelf life

Store dry at <25°C (78°F). Do not use after expiry date.

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