LiWa®

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¿Qué es LiWa?

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Instructions for use

What is LiWa?

LiWa is a ready-to-use, light-curing carving compound of firm, wax-like consistency and consists of acrylates, waxy constituents, fillers and light activators that keep polymerisation contracti-on to an absolute minimum.

How should LiWa be stored?

At a high room temperature (on hot summer days >27°C) LiWa should be stored in a cool and dark place. This will extend the processing life and prevent compound “sticking”. Always keep LiWa in the light-tight original container/packaging. Only remove as much LiWa compound as is actually required. LiWa must not be exposed to natural or artificial light for a prolonged period as this would cause it to polymerise.

How is LiWa cured?

LiWa curing is by light polymerisation. The light activators have been adjusted in such a manner that curing is activated by UV light up to the visible range, i.e. within the range of 280 – 520 nm. Any customary laboratory light curing device, be it a halogen or UV lamp or a combination of the two, can be used for this purpose. Adequate cooling of the object during polymerisation must be ensured.

How much time does it take to cure LiWa?

The curing time depends on the type of lamp, the lamp’s power and age, and the colour of the gypsum model. Dark gypsum reflects less light than a bright gypsum colour. Fluorescent constituents in the gypsum shorten the polymerisation time. A LiWa thickness of 2-3 mm requires an average curing time of 15 – 20 seconds with the Heraeus-Kulzer® Uni XS lamp. (From the fastest light activator Uni XS 30 sec. To the normal instrument approx. 2 – 4 min.).

How is the correct polymerisation time established?

The following test has proved to be successful: Apply an approx. 3 cm long cord, 2-3 mm in diameter, of LiWa on a customary insulated laboratory gypsum model. The base of the cord should be well adapted to the insulated model. After a previously defined polymerisation time (e.g. 3 minutes) check curing by removing the cord approx. 2-3 mm from the model with a thin cut-off wheel. The parted section must be fully cured throughout. If there are soft points inside the core or on the model side, then the polymerisation time has to be lengthened and a new cutting test performed until complete curing has been achieved. Please note: The type and intensity of the light, the duration, colour of the gypsum, LiWa curing time and the age of the lamp all play an important role. Consequently, carry out a test polymerisation regularly in the above described manner within the framework of quality assurance.

How does LiWa behave after curing?

After light polymerisation, LiWa can be mechanically processed as a plastic with rotating tools such as milling, diamond grinding, abrasive paper and rubber polishing. LiWa is sufficiently stable so as not to become deformed, yet sufficiently flexible to lift off a clamp on the underneath area. The recovering ability of cured LiWa is excellent, thus making it one the noteworthy strengths of this new compound.

How is LiWa processed?

LiWa can be “hot” or “cold” carved (i.e. by kneading and working). The easiest way to complete fusing for ho carving is to use an electrical wax knife at a temperature of 150 to 250°С. If the compound starts to emit vapours when the wax knife is inserted in the mas, then progressively diminish the temperature. If the electrical wax knife sticks, then increase the temperature. LiWa behaves like dead wax, i.e. it does not recover, but nor does it contract. Fusion causes LiWa to flow into a thin layer, thereby producing a smooth, homogenous surface that does not perceptibly change after polymerisation.

Since LiWa remains viscous longer it can be cooled down with a cold spray (LiWa-Cool) to produce a consistency that is more readily processed. LiWa can be plastically formed by cold carving (i.e. at room temperature). The consistency is comparable to that of a kneading compound. Working with the fingers achieves the required plastic consistency. Protective gloves must be worn, as prescribed by the professional association, when working with acrylates and materials containing acrylates in order to prevent a possible sensitising or allergic reaction. Wet the fingers with water to prevent compression of the carved sections. This is also beneficial when cooling the carving compound. As is customary in plastics facing, pointed metal and silicone instruments can be very useful in cold carving. Furthermore the basic shape should be slightly over-dimensioned during cold carving, and then finally shaped with rotating instruments after light polymerisation has been completed.

The customary carving waxes are used to add minor supplements and to finish fine carving.

How is gypsum or plastic insulated from LiWa?

Gypsum insulation: LiWa Iso Step I seals the gypsum surface without coating (diffused into the gypsum pores). Action time at least 50 seconds. Then insulate with LiWa Iso Step II (bow away any excess). Action time at least 30 seconds. Plastic insulation: The model is insulated twice in succession with the new LiWa Universal-Isolat. Action time after each application 30 seconds.

LiWa®

Instructions for use

Can also be used to insulate LiWa from LiWa, e.g. sliding parts, wax, plastics and materials.

Important! Model sections/stumps must not be blocked out with a light-curing blocking-out material, nor must it be finished with a light-curing dull varnish. Nevertheless, in spite of careful insulation after LiWa polymerisation, an indetachable connection can be formed with the light-curing-blocking-out material or dull varnish.

How can LiWa be blocked out?

The necessary blocking-out work for the model casting can be completed in the customary manner with tin foil or wax, and this can be precisely adapted to the course of the jaw. It, too, has to be pretreated with LiWa Universal-Isolat 1x.

How can LiWa parts be joined together?

Uncured LiWa parts can be easily joined by fusion with an electric wax knife. Parts already polymerised and ground must be activated at the joint with LiWa Akti Bond.

How can an inhibiting layer (minor smear film) be removed?

Just as any other light-curing material or light-curing material combination, LiWa forms in conjunction with air/oxygen an inhibiting surface layer after light polymerisation. This smear film can be easily removed with alcohol or monomer.

How can a trimmed LiWa surface be dispersion hardened/burnished?

A compact, mirror-smooth surface is produced by applying a thin layer of LiWa Finish-Varnish (a clear surface sealing varnish). This simplifies removal and at the same time trimming time in metal is significantly reduced. The compact layer of LiWa Finish-Varnish gives delicate LiWa work added stability. LiWa finish-Varnish is also ideal for surface sealing of wax models, strengthening contact points and as a retention adhesive.
Important! LiWa Finish-Varnish must only be partially applied on large points/sections that are to be sealed because, otherwise, [with large area surface treatment] light polymerization can delay the entire course of work. The transparency of the varnish limits the polymerisation time to a minimum.

How can LiWa be joined?

LiWa can be joined together with the 2 and 4 mm round wires included with the system. However, conventional wax wires can also be used. The well known, generally valid rules regarding the attachment point of the wires and the casting of the wax objects (depending on the casting unit) must be observed in this context.

How can LiWa be invested?

In principle LiWa is a plastic that is combusted without any residues when preheated. Just as with all plastics, LiWa tends to swell during the heating stage. Consequently, special importance must be attributed to the investment. It must be stable and should not tear when LiWa swells. However, these problems rarely arise with controlled heating, but a few specimen castings should be made with speed casting to find the optimal investment. Since an investment model can be dispensed with, it is not only suitable as a model casting investment but also as conventional K&B investments. These have the advantage of a finer surface. Excellent experiences have also been made with, e.g. Heravest Speed (Heraeus-Kulzer), Formula fast (Dentona AG) and Vulcan (HM Dentalbedarf).

Warranty exclusion:

The fore-going details have been established by numerous laboratory and practical tests. However, the user is responsible for correct application and processing. Warranty claims cannot be made for faults arising in the course of the user’s work.

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