D

Light-Bloc®

Gebrauchsinformation

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Light-Bloc®

Instructions for use

Light-Bloc is a light-curing acrylate-containing plastic for model preparation in dental technology. The plastic carving material is supplied in syringes.

Composition

Light-Bloc Paste consists of multifunction methacrylates with micro- and macro-fillers. Dentaseal contains oligomeric dimethacrylate.

Applications

Light-Bloc material is used for model preparation in dental technology. It is used to compensate intersecting sections, and to build-up and supplement stumps and models.

Processing

Press the paste out of the syringe and remove with customary laboratory instruments and apply onto the model or stump. Light-Bloc Paste is modelled and smoothed with a brush or instrument. If Light-Bloc Paste is applied the plaster should not be moistened. A UV light source should be preferred for polymerisation, but a halogen light source can also be used. After the material has been cured it can be processed with a milling cutter.

Storage

Light-Bloc Paste and Dentaseal have a shelf-life of at least three years if they are kept in their original containers and stored at temperatures not exceeding 22°C. This also applies to packs that have been opened and then re-closed. Avoid exposure to direct sunlight.

Application remarks

The application instructions and other precautionary measures described in the directions for the use of this product must be

observed. If patients still develop undesirable reactions or other effects, even though these instructions have been carefully observed, please notify us in all cases and give a precise description of the accompanying circumstances.

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Willmann & Pein GmbH