I

C-Bond®

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C-Bond®

Instructions for use

C-Bond is an adhesive with outstanding adhesion to dentin and etched enamel. It improves adhesion between the composite and tooth surface.

Indication:

C-Bond is used to improve the adhesion of composites on enamel and dentin. Further, it improves the adhesion between glass-ionomer underfill (lining) cement and composite fillings ("sandwich technique").

Application:

The cavity should be prepared according to the rules for adhesive filling technique. Depending on the situation, apply calcium hydroxide paste (Calcident 450) or an underfill (lining) cement in deep cavities. The enamel is etched, rinsed and then dried in accordance with the etch bonding technique.

The glass-ionomer cement filling should not be etched, only roughened.

Once the prepared cavity is dry, 1-2 drops of light-curing C-Bond are dripped onto the mixing pad (immediately reclose the bottle) and a thin layer is applied to the tooth using a brush or foam pellet. Excess material is removed using a dry brush or a pellet. The material is then cured using halogen light (15 sec.).

Side effects:

Undesirable side-effects are extremely rare if the product is prepared and used properly. Should you become aware of any undesirable side-effects, including in cases of doubt, please notify us, giving as precise a description of the attendant circumstances and symptoms as possible.

Working instructions:

Ensure that the working instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and

give a precise description of the accompanying circumstances.

Contra-indications:

Do not use C-Bond in the area of the pulpa.

Storage/shelf-life:

3 years, store in a cool, dry place. Protect against light!

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Willmann & Pein GmbH