## BH

# Securafil®

### Uputstvo za upotrebu

Securafil je stakloionomerni cementni materijal za punjenje dizajniran za upotrebu u klasi III i V šupljina, kao što su erozije u obliku slova V, karijes u korijenu, karijes na ivicama krunica, karijes na aproksimalnim tačkama kontakta (pripremanje tunela) i u zubnom mesu. Securafil se može također kod zatvaranja zubnih pukotina koristiti. Securafil se veže sa dentinom i caklinom pomoću reakcije karboksilne grupe polikarbonske kiseline sa kalcijem hidroksiapatita i kolagenom. Ova hemijska reakcija čuva mnogo više zubnog materijala tako da bilo kakva podrezivanja zuba pri dnu nisu neophodna. Stakloionomerno punjenje sadži fluorid. Securafil je dostupan u sljedećim nijansama boje: A2, A3

#### Priprema kanala

Zdrav zubni materijal bi trebao biti sačuvan što više da nikakva podrezivanja zuba zbog boljeg prionjanja ne budu potrebna. Rub zuba ne bi trebao biti tanje od 0,5 mm obraden.

#### Zaštita pulpe

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# Securafil®

### Instructions for use

Securafil is a glass inomer cement filling material designed for use in class III and V cavities, such as V-shaped erosions, root caries, caries in crown margins, caries in the approximal contact points (funnel preparation), and in delicious teeth. It may also be used as fissure sealant. Securafil bonds to dentine and enamel by reaction of the carboxic acid with hydroxy apatite calium and with collagen. This chemical reaction preserves sound material. No undercuts are necessary. The glass inomer filler contains fluoride. Securafil is available in follwing colourshades: A2, A3.

#### Cavity preparation

Sound dentine should be preserved as far as possible since undercuts are not necessary due to Securafil's good adhesion. Cavity design preferably should allow the edge thickness of the cement layer to be more than 0,5 mm.

#### Pulp protection

Cavities are rinsed carefully with water and dried with oil free air. Deep cavities should be protected with a calciumhydroxide liner.

#### Application of cleanser

Apply tooth cleanser on freshly prepared dentine surfaces or use it to remove plaque. Leave the tooth cleanser on the surfaces for 10 seconds, rinse with water and dry using oil free air. Avoid Dehydration!

#### Mixing of Securafil

Shake powder before dispensing. Use two level measuring scoops of powder for two full drops of liquid (without air-bubbles, hold bottle upright, press bottle slightly). Mix powder into liquid beginning with half of the powder and continue taking half of the remaining powder and so on. To minimize evaporation of the liquid use only a small area on the mixing pad.

#### Application of Securafil

The ready mixed material is applied into the cavity. The material is condensed slightly. The surface is shaped using a matrix. If used without matrix, it should be covered with an inomer varnish. Securafil is malleable for 2 min (incl. Mixing). Once the hardening process starts, the cement should not be disturbed. Setting time is 3 min.

#### Finishing

After the setting of the matrix can be removed. The now unprotected area should be covered with an inomer varnish. It is recommended to finish the material after additional 6 min. Excess of material can be removed with white stones or flexible discs which have been lubricated with petroleum jelly. Note: Do not use a waterspray!

Cover the whole restoration with inomer varnish to protect the inomer surface from water for 24 hours.

**Side effects** Undiserable side-effects to this medical-product are not to expect, if working corres- ponding to the instructions of use. In case of side-effects - even when doubtful - please inform us with exact description circum- stances and symptoms.

#### Instruction references

Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made according to the Instructions of use) please inform us, if possible, with exact description of the accom- panying circumstances and symptoms.

#### Storage

Stored at room temperature under dry conditions, the shelf life is 3 years if the packing is unopened and 6 months when the package is opened.

#### Precautions

Avoid contact of unset material, varnish and tooth cleaner with the eyes.

**Note:** After the application clean the instruments under water immediately.

**Presentation** [REF]WP4142 - Securafil, A2, 15g Powder / 10g liquid + measuring spoon [REF]WP4143 - Securafil, A3, 15g Powder / 10g liquid + measuring spoon [REF]WP4095 - Securafil, A2, 7g Powder / 5g liquid + measuring spoon [REF]WP4096 - Securafil, A3, 7g Powder / 5g liquid + measuring spoon [REF]WP4145 - Securafil, 10g liquid

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