D

## Ideal Light®

### Gebrauchsinformation

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## Ideal Light®

### Instructions for use

Low viscosity, addition-curing, elastomeric precision impression material based on polyvinyl siloxane – with an especially short intra-oral setting-time.  
Type 4823, Type 3; light-bodied consistency

#### Range of indications

- Syringeable impression material for putty/wash impressions (one step and two step), sandwich impressions and putty/wash impressions taken by using a foil.

- Low viscosity impression material for relinie impressions.

#### Working and setting times

Total working time:	2 Min
Intraoral setting time:	3 Min
Setting time:	5 Min

**To ensure excellent impressions, the total working time and intraoral setting time must be adhered to exactly. Should the temperature of the material deviate noticeable from 23° (74°F), the total working and setting times will be affected.**

#### Instructions for the use with the dispensing gun

- Insert the cartridge into the dispensing gun. Ensure that the notches on the cartridge base are pointing downwards. The clip will not close if the cartridge is not properly inserted into the delivery system.
- 1) Remove the cartridge closure.
  - 2) Do not mount the mixing tip. Extrude impression material until it exits both chambers at the same rate.
  - 3) Mount the mixing tip on the cartridge and extrude the required amount of material, e. g. directly into the preliminary impression, or into a delivery syringe.
  - 4) Leave the filled mixing tip on the cartridge after use or apply the cartridge closure.
  - 5) Before using the cartridge again, remove and discard the used mixing tip/the cartridge closure. Check the cartridge openings to ensure that no plugs of hardened material remain. Should this be the case, unblock the openings with a pointed instrument and check that the material is extruded from the both at the same rate.
  - 6) Mount a fresh mixing tip and proceed as usual.

#### Technical data

Storage conditions: In a dry place, at 18-25°C (64-77°F).  
Disinfection: Impression can be disinfected.  
Electroplating: Such impressions may be copper or silver

plated.  
Pouring the impressions: After 1 hour following disinfection before pouring the impression, the impression can still be poured up to 14 days after time  
Shelf life of impression material: See printing  
Linear dimensional change: -0,20%  
Recovery from deformation: 99,7%  
Strain in compression: 5,0%  
Shore A hardness: 47

#### Warning

Prior to taking the wash impression, the preliminary impression must be cleaned, dried and trimmed, using standard methods.

Solutions applied prior to impression taking (such as retraction solution) may impede the setting reaction of the impression material and must therefore be cleaned off by spraying thoroughly with water and drying.

Do not leave remains of the impression material in the sulcus.  
Avoid splattering this material onto clothing as it cannot be removed by dry cleaning. The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded.

#### Model fabrication

For model fabrication do not pour the impression until 1 hour after disinfection.

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