BH

P-Cem® Tuben

Uputstvo za upotrebu

P-Cem je cement bez eugenola koji je iznimno lagan za oblikovanje i koji ne inhibira stvrdnjavanje umjetnih materijala za punjenje i pričvršćivanje. U slučaju alergija na eugenol preporučljivo je koristiti P-Cem. P-Cem sadrži kalcijum hidroksid i hidroksilapatit za zaštitu i očuvanje zuba.

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Verpackung	Side Effects
REF WP4010 25 g Basis-Paste in Tube	Undiserable side effects of this medical-product are not to expect if working corresponding to the instruction of use. In case of allergies to eugenol – even when doubtful – please ecription of accompanying circumstances and symptoms.
25 g Katalysator-Paste in Tube	Instruction references
	Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made to the instruction of use) please inform us, if possible, with exact description of accompanying circumstances and symptoms.

GB

P-Cem® Tuben

Instructions for use

P-Cem is an eugenol-free easily to process cement. It does not influence the setting of resin-based filling and luting materials. It is recommended to use P-Cem in case of allergies to eugenol. To support the preventies of tooth vitality P-Cem contains calcium hydroxide and hydroxylapatit.

Field of Application	Side Effects
P-Cem is used as a temporary cement for attachment and sealing.	Undiserable side effects of this medical-product are not to expect if working corresponding to the instruction of use. In case of allergies to eugenol – even when doubtful – please ecription of accompanying circumstances and symptoms.
Method of use	Instruction references
P-Cem should be used at room temperature. In case of lower temperature (use of a cooled glass-plate) the working time will be extended correspondingly. Base and catalyst paste should be mixed in proportion: 2 : 1 v/v. For that purpose, put one length of the catalyst paste and accordingly double the length of base on a mixing pad of strong, non-absorbing material or glass-plate respectively (1 : 1 w/w). Then mix the two pastes for approximately 30 seconds until the material takes on a homogeneous consistency. Use a plastic spatula.	Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made to the instruction of use) please inform us, if possible, with exact description of accompanying circumstances and symptoms.

Apply a thin layer inside the temporary crown or bridge and affix to the prepared tooth. Do not use dry cores, moisten them with some water beforehand. P-Cem remains workable for about 1-2 minutes and sets in the mouth within 3-6 minutes. 3-5 minutes after P-Cem has set, remove any excess.

A humid milieu (humid cores) additionally facilitates to remove work afterwards. A harder and faster setting material is obtained by adding more catalyst paste (temporary sealing), or zinc oxide. Extensive superficial and multisurfaced cavities as well as temporary longtem sealing work (max. 1 month) should be treated with a harder cement. Residues of P-Cem can be removed from crowns easily after they have been shortly sprayed with a cooling spray.

Presentation	Side Effects
REF WP4010 25 g base in tube,	Undiserable side effects of this medical-product are not to expect if working corresponding to the instruction of use. In case of allergies to eugenol – even when doubtful – please ecription of accompanying circumstances and symptoms.
25 g catalyst in tube	Instruction references
	Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made to the instruction of use) please inform us, if possible, with exact description of accompanying circumstances and symptoms.

Do not exchange the close-ups of the tubes.

Apply a thin layer inside the temporary crown or bridge and affix to the prepared tooth. Do not use dry cores, moisten them with some water beforehand. P-Cem remains workable for about 1-2 minutes and sets in the mouth within 3-6 minutes. 3-5 minutes after P-Cem has set, remove any excess.

A humid milieu (humid cores) additionally facilitates to remove work afterwards. A harder and faster setting material is obtained by adding more catalyst paste (temporary sealing), or zinc oxide. Extensive superficial and multisurfaced cavities as well as temporary longtem sealing work (max. 1 month) should be treated with a harder cement. Residues of P-Cem can be removed from crowns easily after they have been shortly sprayed with a cooling spray.

Supplementary indications

To obtain a softer material and to keep P-Cem workable for a longer time, add a little bit more of P-Cem-base or a small trace of Vaseline. A humid milieu (humid cores) additionally facilitates to remove work afterwards. A harder and faster setting material is obtained by adding more catalyst paste (temporary sealing), or zinc oxide. Extensive superficial and multisurfaced cavities as well as temporary longtem sealing work (max. 1 month) should be treated with a harder cement. Residues of P-Cem can be removed from crowns easily after they have been shortly sprayed with a cooling spray.

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