I

Calcident LC®

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Calcident LC[®]

Instructions for use

Calcident LC is a light curing radiopaque one-component material containing calciumhydroxide.

Indications

- Indirect pulp capping
- Liner under all filling materials
- Protection against acid and when using the etch technique

Application

The dental surfaces have to be dried thoroughly and contamination avoided. The preparation corresponds to the rules of the utilized filling technique. Application of an aqueous calcium hydroxide preparation is recommended in close proximity to the pulp. To prevent marginal leakage and infiltration of bacteria in the cervical/approximal margins, permanent insulation is recommended. Calcident LC can be directly applied from the syringe. Cure Calcident LC with halogen light. Keep the light source as close as possible to the surface (approx. 1-2mm) and cure for at least 40-60s. Carefully remove excess material with a rotating instrument and apply the filling according to the manufacturer's instructions.

Caution

Light devices should emit at 450 nm and be regularly inspected. The light should be placed as close as possible to the filling material.

Contra-indications

Not for use at patients who are known to be allergic to derivatives of methacrylic acid.

Side effects

No side affects are known to date. Hypersensitivity against components of the material cannot be excluded in every case. Avoid skin contact with the paste. In the event of inadvertent contact, wash immediately with soap and water. In the event of eye contact, immediately rinse with plenty of water and seek medical treatment. If you experience any undesired effects, even in case of doubt, we ask for notification in every case, describing the concomitants and symptoms as detailed as possible.

Processing instructions

Observe the processing instructions in the Instructions for use included with the product and the precautionary measures under all circumstances. Should the patient still show an undesired reaction or if there are undesired effects, we ask for notification in every case, describing the concomitants as detailed as possible.

Warnings/precautions

Calcident LC adheres to materials, e. g. crown and bridge materials, temporary sealants on methacrylate basis, certain impression materials (e.g. polyether). This could lead to the removal of the lining (core build-up). In such cases it is recommended to apply an insulating agent (e.g. vaseline, glycerine gel or similar). Remove insulating agent prior to definite application. The consistency of Calcident LC does not allow to remove excess material with a probe, this could eventually remove the whole lining (use a rotating instrument).

After use Calcident LC avoid direct light (operation lamp, sunlight) and apply the material immediately. Calcident LC contains calcium hydroxide, TEGDMA , BHT, UDMA and amine. Do not use at patients with known allergies to these ingredients. Eugenol and other contaminates (e.g. thymol prevent the polymerization and therefore should be avoided. Insufficient curing times and/or light intensity lead to incomplete curing. Therefore the halogen lamp/light guide should be checked regularly. Tightly close after use and store in an area protected from light and moisture at temperatures between 0°C – 22°C.

Note

Applicator-tools (tips, brushes) are only determinate for single-use at only one patient.

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