WP PerfectPrint -Tray-®

WP PerfectPrint -Tray-®

Gebrauchsinformation

Allgemeine Informationen

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General information

WP PerfectPrint -Tray- is a methacrylate-photoinitiators, pigments and additives, based light-curing resin for use in the dental laboratory. For DLP 3D printers operating at a wavelength of 385-405 nm.

Intended use

Production of dental models for 3D printing technology

Special precautions

- Non-compliance with instructions and improper use may lead to impairment of quality.
- Wear protective gloves, protective clothing and safety goggles while working.
- Avoid contact with eyes! If liquid material comes into contact with eyes inadvertently, rinse immediately with plenty of water and consult a doctor if necessary.
- Avoid direct skin contact with liquid material and components prior to post-curing.
- Keep out of reach of children.
- Observe safety data sheet.

Side effects

- Irritation of eyes and skin.
- The polymerised material should not be used if there are allergies to any of the ingredients or in case of contact allergies.
- In rare cases, allergic reactions to components of WP PerfectPrint -Tray- may occur. Inhalation of monomer vapour and skin contact with non-polymerised materials should be avoided.

Notes on use

- Shake the bottle for at least one minute before use.
- When filling the pan with resin, observe the instructions of the device manufacturer.
- Note the temperature that the device manufacturer specifies for production. The resin can also be processed at room temperature.
- Darken the room slightly before removing the print platform with the printed object. Then, open the hood and remove the printing platform with the printed object.
- Remove the printed object from the platform.
- Only remove the support structures after cleaning and post-curing.

Using used resin

- After sieving (removal of polymerised particles) with a sieve gauze (about 200 µm), the used resin can be mixed with fresh resin 1:1 and reused.

Post-processing of the printed object

- Clean the printed object with ethanol or with a soft brush for a maximum of 3 minutes.
- Optional: Clean printed object in a separate container with fresh ethanol in an ultrasonic bath as needed (5-10 mins).
- Optional: Clean openings, holes and gaps with compressed air.

Post-exposure of the printed object

It is recommended to post-cure the printed object on each side for 20 minutes in a UVA light-curing device with 4 lamps of 9 watts each. Please note the information provided by the device manufacturer.

Restrictions on use

- Intended exclusively for the creation of dental models in the dental sector.
- Must not be brought into direct contact with the patient.

Instructions for Use

Classification

Class I medical device (classification in accordance with Annex VIII, Chapter III, Rule 5 of Regulation 2017/745).

Composition

Mixture of acrylate and methacrylate resins, photoinitiators, pigments and additives.

Storage and durability

- Store at room temperature/0-22 °C in a dark, dry place.
- Only handle WP PerfectPrint -Tray- as briefly as possible to reduce exposure to daylight/room light. Excessive exposure to light leads to an unwanted hardening of the material.
- Close the container immediately after each use and keep it tightly closed.
- Do not use after the expiration date!

Packaging

- WP5189 WP PerfectPrint -Tray-, orange, 1000 g
- WP5190 WP PerfectPrint -Tray-, orange, 500 g
- WP5191 WP PerfectPrint -Tray-, blue, 1000 g
- WP5192 WP PerfectPrint -Tray-, blue, 500 g
- WP5193 WP PerfectPrint -Tray-, pink, 1000 g

Reporting obligation

Any serious incident that occurs in connection with the product, such as the death or temporary or permanent deterioration in the health condition of a patient, user or other person, and any serious risk to public health that occurs or could have occurred in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.