

Extra-Gel

Description:

Extra-Gel is a thixotropic etching gel with 37 % of phosphoric acid. The product comes in a user-friendly packaging in a one-chamber syringe with a very thin tip. Thanks to that, Extra-Gel can be dosed precisely and applied at the aimed position.

Properties:

- 37 % phosphoric acid
- Thixotropic

Advantages:

- Safe and precise dosage possible
- Point-precise application
- No uncontrolled flowing

Indication:

Recommended for etching the tooth surface as part of the total-etch technique for direct and indirect Restorations with Composite.

Physical Data:

Density: 1,3 – 1,4 g/ml

Storage conditions: dry and protected from light at 0 – 22 °C

Shelf life: 3 years

Presentation:

WP4050: Extra Gel, 1 x 2,0 ml syringe + blue cannulas

WP4051: Extra Gel, 2 x 2,0 ml syringe + blue cannulas

