

Fisseal

Description:

Fisseal is a light-cure composite for fissure sealing. The material helps prevent caries in the occlusal surfaces. A particular production technique ensures that the material adjusts very well to the fissures. This guarantees long lasting fissure sealing. Fisseal is available in two variants, transparent and white. The white material additionally strengthens the enamel by continuously releasing fluoride.

Properties:

- Light-cure
- Available in two shades: white and transparent
- The white material contains fluorides

Advantages:

- Time saving thanks to light-curing
- Caries prophylaxis
- Doubled caries protection of the white material thanks to the release of fluorides
- An esthetic solution with the transparent material

Indication:

Suitable for occlusal fissure sealing for a permanent caries prevention.

Physical data:

Colour: white and transparent
Curing: after 3 Min.
Viscosity: 80 - 100 dPa

Storage conditions: dry and protected from light, at 0 – 22 °C.

Shelf life: 3 years

Presentation:

WP4170: Fisseal white, 2 x 2 ml syringes + 1 x 2 ml syringe Extra-Gel + cannulas

WP4171: Fisseal transparent, 2 x 2 ml syringes + 1 x 2 ml syringe Extra-Gel + cannulas

