

Description:

UVA light curing unit / tunnel version for the polymerisation of light-curing tray plates. The unit has a switch with two preset times (120 sec. / 180 sec.).

Features:

- LC LED box
- Display with operating time
- Four preset times (1 min. / 3 min. / 5 min. / 10 min.)

Advantages:

- Easy handling
- Complete curing of tray material

Indication:

For curing the light-curing tray plates Plaque Photo.

Physical data:

Adapter: Rated input power: 100-240 VAC 50/60 Hz.

Rated output power: 24 VDC/ 1.5 A

Maximum consumption: max. 48 W

Lamp: 36 LEDs

Wavelength: 365+405nm

Dimensions (W/D/H): 245 x 200 x 100 mm

Weight: 905g

Presentation:

WP6095: Highlight Easy LED

