

LiWa Universal Isolant

Description:

LiWa Universal Isolant is a liquid that insulates LiWa modellations against other materials, like steel, wax, resins or composites and also to insulate LiWa modellations against each other. LiWa Universal Isolant is self-curing and should be applied in two layers.

LiWa Universal Isolant is part of a large range of resins in various consistencies, preformed resin patterns, useful complementary materials for different areas of use and practical sets.

Properties:

- Liquid
- Self-curing
- Easy to handle
- Application of two layers necessary

Advantages:

- Time-saving
- Good isolation of LiWa against other materials, like resins or steel

Indication:

Liquid for the isolation of LiWa pieces against other materials, like steel or resin.

Storage conditions: dry and protected from light at 0 - 22 °C

Shelf life: 3 years

Presentation:

WP5083: LiWa Universal Isolant, liquid, 10 ml bottle

