

**Description:**

Compact, electric waxer with infinitely variable temperature adjustment for working with LiWa light-curing, burnable resins. LiWaxer® has an ergonomically shaped handpiece with easy-to-change modelling tips. The modelling awl is always included; three additional modelling tips are available as optional accessories (blade, spatula, knife). The heating phase is short and takes only 20 to 30 seconds at 7.5 watts.

**Features:**

- Electric wax knife with modelling awl
- 20 - 30 sec. heating phase
- 4 different modelling tips available

**Advantages:**

- Short heating phase
- Stepless temperature adjustment
- Enables quick working with LiWa

**Indication:**

For "hot" modelling of the light-curing, burnable resins LiWa, LiWa II and LiWa Form.

**Physical data:**

Dimensions (L/W/H): 100 x 65 x 55 mm  
Operating voltage: 220-240 V / 50-60 Hz  
Piston / handpiece: 12 V / 7.5 W  
Stepless regulation: from 50°C - 350°C  
Heating control lamp: LED

**Presentation:**

WP6110: 1 LiWaxer modelling station + modelling awl + plunger tray  
WP5105: Spare tip (awl)  
WP5106: spare tip (knife)  
WP5107: spare tip (blade)  
WP5108: spare tip (spatula)  
WP5109: Spare flask

