

Pro-Dye LC® / Pro-Dye Thinner®

Description:

Pro-Dye LC is a light-curing varnish for cores in the model. The material is a one-component liquid, delivered ready to use in the colours red, blue, green, yellow and nature. After the polymerization with intense UV or halogen light, this core varnish builds a hard and smooth surface. Hence it builds a scratch-proof surface within seconds.

Pro-Dye Thinner dilutes the core varnish without changing its properties after curing.

Properties:

- Liquid consistency
- Light-curing with UV or halogen light
- Scratch-proof and smooth surface
- Core varnish Pro-Dye LC is available in 5 different colours
- Pro-Dye Thinner dilutes Pro-Dye LC, without changing its properties after curing

Advantages:

- Easy handling
- Quick to use
- Hard, scratch-proof and smooth surface
- Opaque colours
- Shiny surface, after removing the smear layer

Indication:

The core varnish Pro-Dye LC is used as spacer between the core in the model and the crown.

Physical properties:

Layer thickness: 10 my

Storage conditions: dry and protected from light at 0 - 22 °C

Shelf life: 3 years

Presentation:

- WP5518: Pro-Dye LC, blue, 30 ml in black brush-bottle
WP5519: Pro-Dye LC, red, 30 ml in black brush-bottle
WP5520: Pro-Dye LC, green, 30 ml in black brush-bottle
WP5521: Pro-Dye LC, yellow, 30 ml in black brush-bottle
WP5522: Pro-Dye LC, nature, 30 ml in black brush-bottle
WP5517: Pro-Dye Thinner, thinner liquid, 10 ml bottle

