

Re-Cord

Description:

Re-Cord is a choice of non-impregnated knitted retraction cords made of 100 % cotton-wool. Re-Cord consists of thousands of interconnecting tiny loops building a tubular cord. It guarantees an easy application, doesn't frizzle and has got a high absorption capacity. Re-Cord remains in the sulcus where it carries-out an optimal pressure. This assures an excellent tissue suppression. Re-Cord is available in five different sizes, separately or in sets of 4.

Properties:

- 100 % knitted cotton-wool
- Non-impregnated
- Even, permanent pressure on the gum groove
- Absorbs large quantities of blood

Advantages:

- Very good retention in the sulcus
- Elastic pressure on the sulcus opening
- Very good hemostatic property
- Does not frizzle

Indication:

Temporary Retraction of marginal Gingiva and therefore presentation of the preparation-border.

Physical data:

Length: 300 cm

Storage conditions: dry and protected from light 0 - 22 °C

Shelf life: 4 years

Presentation:

WP4919: Set with 4 different sizes: size: #00, #0, #1, #2, 4x 300 cm

WP4921: Size: # 000, 1 x 300 cm

WP4915: Size: # 00, 1 x 300 cm

WP4916: Size: # 0, 1 x 300 cm

WP4917: Size: # 1, 1 x 300 cm

WP4918: Size: # 2, 1 x 300 cm

