

# Secura Varnish

## **Description:**

Secura Varnish is an adhesive varnish for sealing the dentin before fitting inlays, gold crowns, jacket crowns or bridge Attachments, and for fixing restoration materials, as well as for the sealing of glass-ionomer-cement and amalgam fillings. Secura Varnish is a clear liquid.

## **Properties:**

- Self-curing
- Liquid texture

## **Advantages:**

- Easy handling
- Esthetic results
- Good sealing of the dentine (e.g. sensitive tooth necks)
- Good sealing of glass ionomer cements

## **Indication:**

Recommended for a straight-forward surface sealing of the dentin (cavities and stumps), of glass ionomer cements (e.g. Securafil) and amalgam fillings, as well as for fixing of sealing materials.

**Storage conditions:** dry and protected from light at 0 – 22 °C

**Shelf life:** 3 years

## **Presentation:**

WP4099: 10 ml bottle

