

# Vision Propy Paste

## **Description:**

Vision Propy Paste is a universal polishing paste that generates excellent polishing results for teeth and acrylics with RDA-value of 237. It is available in two tastes: cinnamon and strawberry.

## **Properties / advantages:**

- Good polishing properties for teeth and composites
- Agreeable tastes of strawberry or cinnamon
- Good bleaching quality for the removal of discolourations

## **Indication:**

- Cleaning before tooth bleaching
- Professional Polishing of the teeth
- Polishing of composites and precious metals
- Polishing of crowns and filling materials

## **Physical data:**

Colour: grey  
Aroma: cinnamon, strawberry  
Solubility in water: 5 g powder / 10 g water  
Density: 1,2 – 1,6

**Storage conditions:** dry and protected from light, at 0-22 °C.

**Shelf life:** 3 years

## **Presentation:**

WP6120: Cinnamon, 1 x 160 g box  
WP6139: Strawberry, 1 x 160 g box



\*Diese Werte basieren auf Literaturrecherche oder externen tests.