

# WP Occlusion

## **Description:**

WP Occlusion is a very fine spray with peppermint odour used to make contact points in prosthetic works visible. The product is only recommended for use in the dental laboratory.

## **Properties:**

- Very fine spray
- Mikro-thin pigmentation film
- Agreeable peppermint odour

## **Advantages:**

- Simple handling
- Easy dosage
- Very good marking of the occlusion
- Easy to remove
- Spray-can can be emptied entirely, as the pigments do not remain in the can

## **Indication:**

To mark the contact points during the making of prosthetic restorations

**Storage conditions:** dry and protected from light at 0 - 22 °C

**Shelf life:** 3 years

## **Presentation:**

WP5110: WP Occlusion, green, spray, 75 ml spray can

