

SMARTER CLEANER FASTER

Pro-Fill®

Light-curing temporary filling material
optimized for clinical efficiency with easy application,
reliable adhesion, and effortless removal



Why Pro-Fill® outperforms traditional solutions:

No Dry-Out, Always Ready

With Pro-Fill, there is no risk of drying out material during storage time. Overall self-curing temporary filling materials do have the problem of drying-out until the material is completely emptied. Thanks to the light-curing properties, Pro-Fill remains stable, is easy to handle and always ready for immediate use.

Application into the cavity can be done quick and straightforward.

Instant Set, No Restrictions

After light activation, the material cures instantly and allows the patient to eat solid food right after treatment – without the risk of the temporary filling coming loose. In contrast, self-curing materials require a longer setting time, which is caused by saliva exposure.

Clean Removal, Faster Workflow

A particularly important advantage of Pro-Fill is that the temporary filling can be removed in one piece by using a suitable instrument. This eliminates the need for laborious cavity cleaning, which is typically necessary when using self-curing temporary filling-materials. As a result, the practitioner can proceed immediately with the next step of the treatment – saving valuable time and increasing overall efficiency.

Packaging

REF WP4041
Pro-Fill, 2 x 4,5 g syringes

Willmann & Pein GmbH

Schusterring 35, 25355 Barmstedt Germany
☎ +49 4123 92280 ✉ info@wp-dental.de

wp-dental.de