



DENT-A-CAV®

We are pleased to inform you that our product Dent-a-cav is now available again – **with improved quality!**



NOW
AVAILABLE
AGAIN

IMPROVED FORMULA.
SAME TRUSTED
RELIABILITY.



Whats new?

- No premature drying when stored correctly
- Cures reliably in the cavity without cracking
- Secure and tight seal that remains intact until further patient treatment

Please note:

Dent-a-cav must always be stored properly (in a dry place at 1°C to 22°C). After dispensing, the container should be immediately closed again and stored under the recommended conditions to maintain product quality over time.





Why Choose Dent-a-cav?


- Non-toxic: gentle on pulp and gingiva
- Excellent marginal seal for temporary cavity fillings
- Resistant to thermal and mechanical stress
- Eugenol-free
- Fast and easy to use



Proven Ingredients. Reliable Results.

Dent-a-cav® is based on calcium sulphate, zinc oxide, and ethyl cellulose.

Its eugenol-free formula makes it fully compatible with all composite restorations — ensuring maximum flexibility for the practitioner.





Simple. Clean. Effective.

Dent-a-cav® is ready to use, applies easily to moist cavities, and hardens in minutes.

No need to fill deep cavities completely — simplifying later removal.

Can be easily removed with a short, strong probe



Available in Two Shades

Choose the right variant for your
workflow or patient preference:

● Pink



REF WP5627 1 × 30 g jar
REF WP4311 2 × 30 g jars

● White



REF WP5628 1 × 30 g jar
REF WP4312 2 × 30 g jars




Need More Information?

If you have any questions about Dent-a-cav®, its use, or availability — don't hesitate to reach out.

We're happy to assist you.

info@wp-dental.de

+49 4123 92280





Willmann & Pein GmbH

Schusterring 35
25355 Barmstedt
Germany

wp-dental.de

Office hours:
Mon. - Thurs. 07:00-16:00
Fri. 07:00-14:30

**Dental Products
and Supplies**