

Calcident LC

Description:

Calcident LC is a light-cure, radiopaque calcium hydroxide paste for indirect pulp capping. Calcident LC can be used as a liner under any filling material. Light-curing takes 40 sec. Calcident LC is available in the one-chamber syringe for the easy application.

Properties:

- Contains 4,2 % calcium hydroxide
- Radiopaque
- Light-curing

Advantages:

- Gentle to the pulp, if it is not exposed and does not show-through
- Fast to use thanks to the syringe and the light-curing

Indication:

Suitable for indirect pulp capping. In direct proximity to the pulp, it is recommendable to use a water based calcium hydroxide material, like Calcident 450.

Physical Data:

Density: 1,1 – 1,2 (g/ml)
Colour: tooth-coloured, opaque

Storage conditions: dry and protected from light, at 0 – 22 °C.

Shelf life: 1,5 years

Presentation:

WP4271: 1 x 2 ml syringe + cannulas
WP4270: 2 x 2 ml syringes + cannulas

