

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

1. Chemical Product and Company Identification

Product name:

Glass Liner

Application of the product:

Light curing lining material

Manufacturer:

Willmann & Pein GmbH

Street:

Schusterring 35

Zip-Code / Town

D-25355 Barmstedt

Contact for technical Information

+49 (0) 4123 9228 0

Phone / Fax / E-Mail

+49 (0) 4123 9228 0 / +49 (0) 4123 9228 49 / info@wp-dental.de

Emergency Information:

Dr. Wolfgang Willmann

Phone: +49 (0) 170 – 405 30 52

2. Hazard Identification

2.1 Classification of the substance or mixture

According to Regulation (EC) No. 1272/2008:

2.2 Label elements

Labelling according to Regulation 1272/2008 (GHS/CLP)

H319 (Causes serious eye irritation)H335 (May cause respiratory irritation)H315 (Causes skin irritation)

Additional hazard advices for human and environment

None

Material Safety Data Sheet

According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

2.2 Label elements

Labelling according to Regulation 1272/2008 (GHS/CLP)



H315

P264
P280
P302+P352
P321
P332+P313
P362+P364

causes skin irritation

Wash...thoroughly after handling.
Wear protective gloves/Protective clothing/eye protection/face protection.
IF ON SKIN: Wash with plenty of water/...
Specific treatment (see...on this label)
If skin irritation occurs: get medical advice/attention
Take off contaminated clothing and wash it before reuse.

H319

P305+P351+P338

P337+P313

causes serious eye irritation

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists: Get medical advice/attention.

H335

P271
P304+P340
P312
P403+P233
P405

may cause respiratory irritation

Use only outdoors or in a well-ventilated area.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.
Call a POISON CENTER/doctor/.../if you feel unwell.
Store in a well-ventilated place. Keep container tightly closed.
Store locked up.

3. Composition/Information on Ingredients

Chemical characterization:

Hazardous ingredients:

Stoff	Cas-Nr./EINECS-Nr./EG-Nr.	Gehalt (%)	Kennzeichnung
Glasionermpulver		0-40 %	H319, H335, H315
1,6-Hexandioldimethacrylate	CAS-Nr. 6606-59-3 EINECS-Nr. 229-551-7	0-10 %	H319, H335, H315, H411
Bisphenol-A-bis (hydroxypropylmethacrylat), Isomere	CAS-Nr. 1565-94-2	0-15 %	H319, H335
4-tert.-Butyl-N,N-dimethylaniline	CAS-Nr. 2909-79-7 EG-Nr. 220-824-6	0-0,10 %	H302, H312, H332, H319, H335, H315
Campherchinon		0-0,10 %	H319, H335, H315

4. First-Aid-Measures

General Advices

Take off with product contaminated clothes.

After skin contact

Wash off skin with soap and water.

Contact a physician if the skin is irritated.

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

After eye contact

Rinse eyes an efficient time.
Contact ophthalmologist if irritation is persisting.

After inhalation

Take person to fresh air.
In the event of symptoms refer for medical treatment.

After swallowing

Seek medical advice immediately.5. **Fire-fighting-measures**

Suitable extinguishing agents

Foam, extinguishing powder, carbon dioxide, water spray.
Adjust fire-fighting-measures to ambiance.

For safety reasons unsuitable extinguishing media

None known

**Particular danger through the product itself,
its combustion products or gases**

None

Protective equipment for fire-fighters

Use conventional measures to combat chemical fires, that means use self-contained breathing apparatus (isolating-unit).

Further information

None

6. Accidental release measures

Personal precautions

Arrange for fresh air.

Environmental protection measures

Do not drain into canalization.

Methods for cleaning / absorption

Take up mechanically.
Apply to local services.

7. Handling and storage

Handling

The conventional precautions for works with chemicals are to be adhered.
Wear overall while working.

Hygiene measures

Keep away from food, drink and animal feedingstuffs. When using do not eat, drink or smoke. Take off contaminated clothing immediately. Avoid contact with skin and eyes.

Information about fire- and explosion protection

None

Storage conditions

Store coll, dry, protected against light and tightly closed.
Store in a well-ventilated place.
Storage time: 3 years if storage conditions will be kept

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

8. Exposure Controls / Personal Protection**Limit values for the exposition at work station and / or biological limit values****Work station limit values (AGW) Germany**

PNCE and DNCL-Values as to be possibly worked out in a preparation safety report do not exist.

Personal protective equipment

Pay attention of protection- and hygienic measures. Protection clothes must be stored separately. The usual precaution measures must be attended.

Breathing protection: not necessary

Eye protection: Wear protection goggles

Hand protection: Sensibilisation: Wearing of applicable protection gloves is recommended. (by experience of the past / literature) for example Ethiparat®-gloves for 5-10 minutes, application of No-touch-Technik, minimal-contact-techniques is recommended

Body protection: Wear business usual protection clothes.

Other protection measures: none

Advices for hygienic

at work station:

Comply business usual protection- and hygienic measures.
Preventive skin protection. Wash hands and skin after finishing work.**9. Physical and Chemical Properties****Appearance**

Appearance: pasty material

Colour: yellowish

Odour/Flavour: acrylates / methacrylates

Safety relevant dataFlash point: not determined
(approx. > 100°C appreciated by the flash point of the ingredients)

PH-value: not determined

Vapor pressure : not determined

Density: not determined

Solubility in water: insoluble

10. Stability and Reactivity**Conditions to avoid**

Protect from heat, sunlight and influence of light.

Cures with exposure to UV light.

Substances to avoid

Free radicals, reducing agents, heavy-metal ions.

Hazardous decomposition products

In case of fire acrid acrylic fumes.

Further information

None

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

11. Toxicological Information**Acute Toxicity**

At skin: irritation
At eyes: irritation
At respiratory system: irritation

Kinetic of the toxic reaction

According the new Material Safety Data Sheet according to REACH, specifications about kinetic of the toxic reaction, metabolism of the substance and allocation of the substance in the human body must be mentioned at this place. Data for this is not available .

Further information

By appropriate exposure and conventionally use and to our experiences and informations the product causes no hazardous effects.

General information

none

12. Ecological Information**General information**

Fish toxicity: not applicable
Water organism: not applicable
Bacteria toxicity: not applicable

Persistence and degradability

Duration: not applicable
Test method: not applicable
Evaluation: not applicable

Further detail to ecology

Do not allow product to reach ground water, body of waters or sewage system.

Bioaccumulation potential

No data available

Effect of the determination of the PTB- and vPvB-Report

Endocrine disruption potential

If the material causes hormonal effects or prevents them is unknown to us.

13. Disposal Considerations**Substance / Preparation**

To dispose of accordingly with notice of the local regulations of authority.
Do not empty into drains.

Disposal key according to AVV regulation

16 03 05 organic wastes containing dangerous substances

Package

Packages that are not contaminated can be recycled.

14. Transport information / Rules**Land transport ADR**

Class: Not a dangerous good according to the transport regulations.
UN-Number: -
Classification code: -
Description of the good: Glass Liner
Packaging group: -
Hazard sheet: -

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

Symbol: -

Rail transport RID

Class: Not a dangerous good according to the transport regulations.
 Classification code: -
 UN-Number: -
 Packaging group: -
 Hazard sheet: -
 Symbol: -

Inland water transport ADN

Class: Not a dangerous good according to the transport regulations.
 Classification code: -
 UN-Number: -
 Packaging group: -
 Hazard sheet: -
 Symbol: -

Maritime transport IMDG

Class: Not a dangerous good according to the transport regulations.
 EmS: -
 Classification code: -
 UN-Number: -
 Packaging group: -
 Hazard sheet: -
 Symbol: -

Air transport IATA

Class: Not a dangerous good according to the transport regulations.
 Packaging-Instruction (cargo) -
 Packaging-Instruction (passenger) -
 UN-Number: -
 Packaging group: -
 Hazard sheet: -
 Symbol: -

15. Regulatory Information**2.1 Classification of the substance or mixture****According to Regulation (EC) No. 1272/2008:****2.2 Label elements****Labelling according to Regulation 1272/2008 (GHS/CLP)**

H319 (Causes serious eye irritation)H335 (May cause respiratory irritation)H315 (Causes skin irritation)

Additional hazard advices for human and environment

None

2.2 Label elements

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

Labelling according to Regulation 1272/2008 (GHS/CLP)



Warning

H315

P264

P280

P302+P352

P321

P332+P313

P362+P364

causes skin irritation

Wash...thoroughly after handling.

Wear protective gloves/Protective clothing/eye protection/face protection.

IF ON SKIN: Wash with plenty of water/...

Specific treatment (see...on this label)

If skin irritation occurs: get medical advice/attention

Take off contaminated clothing and wash it before reuse.

H319

P305+P351+P338

P337+P313

causes serious eye irritation

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists: Get medical advice/attention.

H335

P271

P304+P340

P312

P403+P233

P405

may cause respiratory irritation

Use only outdoors or in a well-ventilated area.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Call a POISON CENTER/doctor/...if you feel unwell.

Store in a well-ventilated place. Keep container tightly closed.

Store locked up.

Hazard decisive components

Contains: -

Safety Assessments

Safety Assessments for substances in this preparation were not executed

National Rules

Water Hazard Class 2 (Estimation on the basis of the ingredients)

Material Safety Data Sheet
According to decree (EC) No. 1907 / 2006

Rev. 4

Revision Date 22.08.2017 SW

16. Further Information**Further applicable EG Guidelines**

Guideline for preparation (1999/45/EC), at last changed by the guideline 2006/8/EC.
Substance Guideline (67/548/EWG), at last changed by the guideline 2006/121 EC.
REACH-Guideline (EC) No. 1907/2006

Manufacturer's suggested application restriction:

For dental use only

R-phrases on which is referred in chapter 2 and 3

R36/37/38 Irritating to eyes, respiratory system and skin

Changes compared to the last version

Change of GHS-Labeling
Change of Composition/Information on Ingredients
Change of Fire-fighting-measures
Change of Handling and storage
Change of Exposure Controls / Personal Protection
Change of Physical and Chemical Properties
Change of Toxicological Information
Change of Ecological Information
Change of Disposal Considerations

Department issuing MSDS:

Dr. Wolfgang Willmann (Phone: +49 (0) 41 23 – 9228 0)

Additional Information:Origin: ¹ <http://www.baua.de>² <http://www.arbeitssicherheit.de>

Some of the information presented and conclusions drawn herein are from sources other than direct test-data on the product itself.

The contents and format of this MSDS are in accordance with EEC Commission (EC) No. 1907 / 2006.

Disclaimer of Liability

The information in this MSDS was obtained from sources which we believe are reliable. However, the information is provided without any warranty, express or implied, regarding its correctness.

The conditions or methods of handling, storage, use or disposal of the product are beyond our control and may be beyond our knowledge.

For this and other reasons, we do not assume responsibility and expressly disclaim liability for loss, damage or expense arising out of or in any way connected with the handling, storage, use or disposal of the product.

This MSDS was prepared and is to be used only for this product. If the product is used as a component in another product, this MSDS information may not be applicable.