

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

1. Chemical product and company identification

Product name:

Securafil LC, Powder

Application of the product:

Glas ionomer filling material, light curing for the dental practise

Manufacturer:

Willmann & Pein GmbH

Street:

Schusterring 35

Zip-Code / Town

D-25355 Barmstedt

Contact for technical information

+49 (0) 41 23 / 9228 0

Phone / Fax / E-Mail

+49 (0) 4123 / 9228 0 / +49 (0) 4123 / 9228 49 / info@wp-dental.de

Emergency Information (24 hours)

Dr. Wolfgang Willmann

Phone: +49 (0) 170 – 405 30 52

2. Hazard Identification

Potential Hazards:

Irritant

Categorization: Xi R36/37/38

Additional hazard advices for human and environment:

none

3. Composition/Information on Ingredients

Chemical characterization:

Preparation of dental glass, benzoylperoxide, pigments.

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

4. First-Aid-Measures

General Advices:

Take off with product contaminated clothes.

Contact a physician if symptoms are obviously showing influence of the product to skin, eyes or through inhalation of dust.

After skin contact:

Wash off skin with soap and water.

Contact a physician if the skin is irritated.

After eye contact:

Rinse eyes an efficient time.

Contact ophthalmologist if irritation is persisting.

After inhalation:

In case of prolonged exposure take person to fresh air.

Seek medical treatment.

After swallowing:

Seek medical advice immediately.

5. Fire-Fighting-Measures

Suitable extinguishing agents:

Waterspray, sand, foam, carbon dioxide.

For safety reasons unsuitable extinguishing media:

Water with full jet.

Particular danger through the product itself, its combustion products or gases:

no information

Protective equipment for fire-fighters:

Use conventional measures to combat chemical fires, that means use self-contained breathing apparatus. (isolating-unit)

6. Accidental release measures

Personal Precautions:

Arrange for fresh air. (well-ventilated areas)

Environmental protection measures:

Do not empty into drains.

Methods for cleaning / absorption:

Take up mechanically.

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

7. Handling and storage

Handling:

The conventional precautions for works with chemicals are to be adhered.
Always keep the bottle tightly closed. Provide adequate room ventilation.

Information about fire- and explosion protection:

none

Storage conditions:

Store material in a tightly closed container at a dark, cool place.
Recommended packaging material: container of glass, metal or plastics.

8. Exposure Controls/Personal Protection

Limit values for the exposition at work station and / or biological limit values work station limit values (AGW) Germany

PNCE and DNCL-Values as to be possibly worked out in a preparation safety report do not exist.

Personal protective equipment:

Personal attention of protection- and hygienic measures. Protection clothes must be stored separately. The usual precaution measures must be attended.

Breathing protection: Do not breathe dust. Dust mask type P2

Eye protection: Protection goggles

Hand protection: Rubber gloves

Body protection: Wear business usual protection clothes

Other protection measures: none

Advice for hygienic
at work station: Comply business usual protection- and hygienic measures.
Preventive skin protection. Wash hands and skin after
finishing work. Do not eat or drink when using.

9. Physical and chemical Properties

Appearance:

Form: powder

Colour: whitish

Odour: weak

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

Safety relevant data:

Melting temperature:	not applicable
Ignition temperature:	not applicable
Flashpoint:	not applicable
Explosion properties:	not applicable
Vapour pressure:	non volatile
Rel. Vapour density:	>1 (methacrylates)
Bulk Density:	0,6g/cm ³
Solubility in water:	insoluble
pH-value:	neutral

10. Stability and Reactivity**Conditions to avoid:**

none

Substances to avoid:

none

Hazardous decomposition products:

Thermal decomposition at ca. 60°C.

No decomposition, if handled and stored according to directions.

Further information:

After considerably exceedance of storage time and / or exceedance of storage temperature it is possible that the product does not have the usual product properties.

11. Toxicological information**Acute toxicity:**

At skin:	irritation
At eyes:	irritation
At respiratory system:	irritation

Kinetic of the toxic reaction:

According the new Material Safety Data Sheet according to REACH, specifications about kinetic of the toxic reaction, metabolism of the substance and allocation of the substance in the human body must be mentioned at this place. Data for this is not available.

Further information:

By appropriate exposure, conventionally use and to our experiences and information the product causes no hazardous effects.

Further details:

May cause sensitization or allergic reactions by skin contact or inhalation

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

12. Ecological information

Fish toxicity: not applicable
Water organism: not applicable
Bacteria toxicity: not applicable
Acute oral toxicity: LD50 (rat) >2000mg/kg

Further detail to ecology:

Do not allow product to reach ground water, body of waters or sewage system!

Bioaccumulation potential:

No data available

Effect of the determination of the PTB- and vPvB-Report:

Endocrine disruption potential

If the material causes hormonal effects or prevents them is unknown to us.

Persistence and degradability:

No data available for the preparation.

13. Disposal considerations

Substance / Preparation

To dispose of accordingly with notice of the local regulations of authority.

Do not empty into drains.

Package

Packages that are not contaminated can be recycled.

14. Transport information/rules

Land transport ADR/RID

Classification (no dangerous cargo)

Class:

UN-Number:

Classification code:

Description of the good: Glass Liner II, powder

Hazard cause:

Packaging:

Packaging group:

Hazard sheet:

Limited Quantities:

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

Maritime transport IMDG/GGVSea

Classification (no dangerous cargo)
IMDG-Code:
UN-Number:
Ems:
Marine Pollutant:
Description of the good: Glass Liner II, powder
Hazard cause:
Packaging
Packaging group:
Hazard sheet:

Air transport ICAO-TI and IATA-DGR

Classification (no dangerous cargo)
Class:
UN-Number:
Description of the good: Glass Liner II, powder
Hazard cause:
Packaging
Packaging group:
Hazard sheet:

15. Regulary Information

Labelling according EU guidelines

Code letter and hazard description of the product:



Xi; Irritant
R36/37/38 Irritating to eyes, respiratory system and skin

Hazard decisive components for the labelling:

Contains: --

National Rules

Water Hazard Class: 0 (own classification)

16. Other information

Further applicable EG Guidelines

Guideline for preparation (1999/45/EC), at last changed by the guideline 2006/8/EC.
Substance Guideline (67/548/EWG), at last changed by the guideline 2006/121 EC.
REACH-Guideline (EC) No. 1907/2006

Manufacturer's suggested application restriction:

Only for commercial use – no public product

R-phrases on which is referred in chapter 2 and 3

36/37/38 Irritating to eyes, respiratory system and skin

Material Safety Data Sheet
according to decree (EC) No. 1907 / 2006

Revision: 3
Revision Date: 02.04.2012 LK

Changes compared to the last version

Change of Hazard identification
Change of Exposure controls/personal protection
Change of Toxicological information
Change of Transport information/rules
Change of Regulatory information

Department issuing MSDS:

Dr. Wolfgang Willmann (Phone: +49 (0) 41 23 – 9228 0)

Additional Information:

Origin: ¹ <http://www.baua.de>
² <http://www.arbeitssicherheit.de>

Some of the information presented and conclusions drawn herein are from sources other than direct test-data on the product itself.
The contents and format of this MSDS are in accordance with EEC Commission (EC) No. 1907 / 2006.

Disclaimer of Liability

The information in this MSDS was obtained from sources which we believe are reliable. However, the information is provided without any warranty, express or implied, regarding its correctness.
The conditions or methods of handling, storage, use or disposal of the product are beyond our control and may be beyond our knowledge.
For this and other reasons, we do not assume responsibility and expressly disclaim liability for loss, damage or expense arising out of or in any way connected with the handling, storage, use or disposal of the product.
This MSDS was prepared and is to be used only for this product. If the product is used as a component in another product, this MSDS information may not be applicable.