

Description:

Securafix is a glass ionomer cement for the permanent fixing of crowns, bridges, inlays, onlays, posts and orthodontic bands. Securafix adheres well to dentin and releases fluoride into the tooth. Securafix is self-curing and radiopaque. It is available as a flexible powder/liquid system.

Properties:

- Self-curing
- Radiopaque
- Contains fluoride

Advantages:

- Fluoride release
- Safe adhesion to enamel and dentin

Indication:

Suitable for the permanent fixing of crowns, bridges, inlays, onlays, posts and orthodontic bands.

Storage conditions: dry and protected from light at 0 - 22 °C.

Physical data: 100 MPa

Shelf life: 3 years

Presentation:

WP4110: 24 g powder, 10 g liquid, dosage spoon

WP4093: 7 g powder, 5 g liquid, dosage spoon

WP4146: 10 g liquid

