

LiWa Cool

Description:

LiWa Cool is a cooling spray which quickly cools LiWa modellations done “hot” with an electrical wax-knife. LiWa Cool is free of CFC and therefore eco-friendly.

LiWa Cool is part of a wide range of materials: resins in different consistencies, preformed resin patterns, useful complementary materials for different areas of use as well as functional sets.

Properties:

- Easy to handle
- Cools quickly down LiWa modellations which have been done „hot“

Advantages:

- Time-saving
- Facilitates working with LiWa

Indication:

Developed to speeds up the cooling of LiWa modellations done with an electrical wax-knife.

Storage conditions: dry and protected from light at 0 - 22 °C

Shelf life: 3 years

Presentation:

WP5079: LiWa Cool, spray, 200 ml spray can

