

Dent-a-cav

Description:

Dent-a-cav is a temporary close-up material, which is delivered as a ready to use paste, in white or pink. Dent-a-cav adheres to enamel and dentin, shows a very good marginal fit and can be easily removed with a probe. Dent-a-cav is gentle to the pulp and the gingiva. Also, thanks to its resistance against chewing pressure, Dent-a-cav allows for a safe close-up. And, as it contains no eugenol, Dent-a-cav does not impede subsequent composite fillings. The material cures quick hardening inside the cavity.

Properties:

- Good adhesion to enamel and dentin
- Safe marginal seal
- Fast curing
- Easy removal
- Containing zinc oxide
- Free from eugenol

Advantages:

- Easy to use
- Fast curing in the cavity
- Good adhesion and excellent marginal seal
- Time saving
- Gentle to the gingiva and the pulp
- Not impeding subsequent composite filling
- Anti-bacterial effect

Indication:

Suitable for temporary close-ups.

Physical data:

Spreadability: 16 -18 mm
Appearance: homogenous paste, white or pink

Storage conditions: dry and protected from light, at 0-22 °C

Shelf life: 2 years

Presentation:

WP4311: 2 x 30 g glasses, pink
WP4312: 2 x 30 g glasses, white
WP5627: 1 x 30 g glasses, pink
WP5628: 1 x 30 g glasses, white

